

smartpft[®]

devices for smart pulmonary function testing



Medical Equipment Europe GmbH

www.meegmbh.com

smartSOFTmee

advanced graphical user interface
based on WIN7 / WIN8 / WIN10 operating systems


database server



interfaces



software architecture



The SQL database, a professional network server solution represents an international standard for data storage.

It guarantees maximum safety for testing and patient data.

Nearly all requirements of data sharing in large networks can be realized.

- ASCII, BDT / GDT
- HL7 *
- PDF (labeled) *
- VPN networking*
- secure eMailing*
- e-Doc printer *
- Hard/Software for BGA devices*
- (* options)

Our modern graphical user interface makes working with the system simple and efficient.

The workflow is user optimized.


Print templates and predicted value equations can be created and modified by very comfortable built-in editors.



smartVsensor®

electronics

special features



This pneumotachograph sensor type is used since more than 35 years for pulmonary function and pulmonary exercise testing because of excellent characteristics.

Low resistance, low dead space, high linearity and absolutely no sensitivity to humidity make this sensor an optimal choice.

The data processing unit is based on a new low-power μ -processor technology.

Digital pressure sensors reduce the risk electrical interferences.

The USB 2.0 PC interface guarantees long-term compatibility for PC system.

The standard software includes the newest [GLI 2012](#) predicted value equations and [smartZeroing®](#).

An online BTPS data sampling module can be added to the USB spirometer.

In all other systems BTPS sensors are integrated as a standard.

[smartZeroing®](#) means, the annoying zeroing happens discreetly and continuously in the background.

smart pft[®] USB

high precision PC based spirometer



PC Spirometer, USB powered


A powerful spirometer with high precision technology and the possibility for upgrading with a wide range of testing options.

The system requires only one USB cable to connect to a Microsoft[®] WIN7, WIN8 or WIN10 based computer.

This spirometer has all connectivity options of an advanced pulmonary function system like ASCII, HL7 and networking.

Standard features are

- Slow vital capacity test
- Flow / Volume test
- PRE/POST comparison
- Trend reporting
- Programmable challenge testing protocols
- [smartZeroing[®]](#)
- Incentive graph for kids
- [GLI 2012](#) predicted values
- [SQL](#) data base



Touch-PC for outdoor applications

This tablet computer device is designed for safe in- and outdoor use.

A massive aluminum housing protects the built in touch computer and other electronic components safely.

This device runs with full testing functionality on batteries even in combination with a shutter system for Rocc and respiratory drive measurements.

All BTPS-sensors are built in.

Options available for all system

MIP, MEP, P0.1 and Rocc
(run without additional power supply)

Rhinomanometry

Compliance test

Table stand

A4 printer

Thermo printer

ASCII-, PDF- and HL7-Interface.

BGA serial and USB interface

smartpft[®] body

massive aluminum & glass fibre construction

combined with

modern design & functionality



Panorama view box

The aim was to create a device that fully meets the requirements of patients.

What makes more sense than constructing an oval shaped box with round glazing so the patient never feels captured.

Very modern and long-term stable materials make it a very useful and unique device.

Technical standard features


Height adjustable, comfortable chair.

3D adjustable device stand. (optionally, e-motor driven)

Ambient pressure, temperature and humidity sensors allow continuous BTPS data sampling and calculation.

The electromagnetic door locking is maintenance free and secure for patients.

The box has a very efficient ambient pressure movement compensation unit.



Designed for patients with poor mobility

The two doors opening concept and a less than 10 cm door step height in combination with handrails facilitate entering this bodyplethysmograph.

Optionally, a special chair mounted on a telescope arm can be added to move wheelchair patients safely into the box.

Also optionally, a 180 kg chair spring is available.

We use only materials with long-term stability

The door rubber is protected by a steel frame.


Box and door frames are made from powder coated aluminum.

Bottom and top plates are made from glass-fiber.

For simple service access, nearly all electronic components are placed on the top and protected by a hood.

smartpft[®] CO-transfer

designed for simplified CO transfer capacity screening



CO transfer capacity test

This device is designed for quick screening of CO-diffusion capacity.

The online gas sampling method allows testing of patients with very low vital capacity starting from approximately 0.7 liters.

The test runs fully automatic and gives online all necessary advices how the patient must cooperate.

The system registers and displays real-time all important test signal data like gas concentrations, volume and mouth pressure.

Standard features are

Selective and fast infrared gas analyzer measuring transfer and inert gas (CO & CH4).

Particular and precise gas supply by a demand valve.

Input mask for Hb values and correlating recalculation of test results.

Automatic altitude corrections by built in ambient pressure sensor.

Online display of mouth pressure during breath holding.

Gas analyzer works without absorber products.

smartpft[®] nebulizer

standardized dosimeter unit for challenge testing



Compressor unit

A comfortable, safe and very robust device to support challenge testing.

No need for an external source of compressed air, a low noise compressor is built in.

The compressor runs fully automatic.

A built in panel computer with LCD-display guarantees a very simple control handling.

Optionally a USB interface allows connection to the PC system.





DeVilbiss 646 nebulizer

Standardized nebulizer type for challenge testing used in many studies.

Using the Rosenthal Chai protocol nebulization method .

This nebulizer and nebulization method fulfill strictly the recommendations of ATS and ERS.

A flow based trigger sensor synchronizes the start of nebulization with beginning of inhalation.

| | | | | | |
|----------------------|--|--------------------------------|---|---|---|
| smartVsensor | materials | housing | PVC | color | white & blue |
| | sensor type | pneumotachograph | screen | variable orifice | hostaphan |
| | flow range | 0,02l/sec - 20 l/s | volume range | 0,02 L to 20 L | numeric integration |
| | dimensions | length = 80 mm | width = 37 mm | height = 25 mm | |
| | backpressure | < 0,055 kPa/l*s @ 15 l/s | linearity error | < 3% absolute | |
| | disinfection: | cold gas & cold liquid | | | |
| pressure transducers | flow | interface I ² C bus | ADC 14 bit | temp-compensated | linearity error < 0,1 %FSS |
| | | type | piezoresistive | differential | |
| | | range | 12,7 | mbar | |
| | mouth pressure | interface I ² C bus | ADC 14 bit | temp-compensated | linearity error < 0,1 %FSS |
| | | type | piezoresistive | differential | |
| | | range | 200 | mbar | |
| | box pressure | interface I ² C bus | ADC 14 bit | temp-compensated | linearity error < 0,1 %FSS |
| | | type | piezoresistive | differential | |
| | | range | 2,5 | mbar | |
| | ambient pressure | interface I ² C bus | ADC 8 bit | temp-compensated | linearity error ± 10 mbar |
| | | type | piezoresistive | absolute | |
| | | range | 500 - 1150 | mbar | |
| | temperature sensor | interface I ² C bus | ADC 8 bit | temp-compensated | linearity error < 1%FSS |
| | | type | PTA internal | range | -40 to + 125 °C |
| | humidity sensor | interface I ² C bus | ADC 8 bit | temp-compensated | linearity error ± 1 %rF |
| | | type | capacitive Polymer | range | 0 - 100 % rF |
| shutter | materials | housing | PVC | color | white & red |
| | type | electromagnetic | | | |
| | effective dead space | < 20 ml | occlusion time | < 60 ms | |
| | dimensions | length = 145 mm | width = 85 mm | height = 85 mm | weight = 450 g |
| | backpressure | < 0,03 kPa/l*s | | | |
| | disinfection: | cold gas & cold liquid | | | |
| smartpft USB | dimensions | length = 163 mm | width = 65 mm | height = 40 mm | weight = 151 gr. |
| | | material | PVC | color | white |
| | | disinfection: | cold gas & cold liquid | | |
| smartpft body | weight | 139 kg with acryl | | | |
| | power supply | 230 VAC | max 350 VA | optional 1000 VA | |
| | materials | door & box frame | aluminum | powder coated | color = white |
| | | base & cover plate | glass fibre | | color = white |
| | | windows | glass / acryl | floor | PVC antrazit |
| | volume | 860 liter | | | |
| | dimensions | with closed doors | length = 115 cm width = 83 cm height = 185 cm |  |  |
| | | with opened doors | length = 120 cm width = 111 cm | | |
| | | | height = 190 cm | | |
| | | | width = 78 cm | | |
| | minimum door dimensions for transport | | | | |
| | box & mouth pressure calibration syringe | | | | |
| | | volume | 50 ml | | |
| | | type | sinus electrical | | |
| | | calibrates | box pressure | mouth pressure | |
| | | voltage | 12 V | | |
| | | current | 700 mA | | |
| | box ventilation valve | | | | |
| | | type | electromagnetic | | |
| | | voltage | 12 V | | |
| | | maximum current | 1,1 A | | |
| | loudspeaker | inside body box | | | |
| | chair | fixed, hight adjustable | | | |
| | | maximum weight | 120 kg | optional | 180 kg |
| | | color | black | material | soft plastic |

technical data II

| | | | | | |
|----------------------|---|------------------|--------------------|-----------------|-------------------------|
| smartpft CO transfer | non dispersive infrared absorption (NDIR) | gas channels | CH4 & CO | | |
| gas analyzer | CH4 & CO channel | range | 0 ppm - 3000 ppm | | |
| | | accuracy | < ± 1% FSS | resolution | < 0,5% FSS |
| | | resolution | < 0,5% FSS | linearity | < 1% FSS |
| | | warm up time | max. 5 min | zero drift | < 1% FSS / 24h |
| tablet touch PC | dimensions | length = 290 mm | width = 171 mm | height = 20 mm | weight = 1350 g |
| | | materials | aluminum | color | RAL 9006 |
| | power supply | voltage | 12 V= | maximum current | 2 A |
| | PC system | processor | Intel® Atom™ Z2760 | screen | 8" WXGA multi-touch LCD |
| | | memory | 2GB | storage | eMMC 32 GB |
| | | interfaces | WiFi | USB | |
| | | operating system | windows 8 | 32 bit | |
| smartpft nebulizer | power supply | | 230V AC / 50 Hz | | |
| | dimensions | length = 250 mm | width = 250 mm | height = 145 mm | weight = 2900 g |
| | reservoir | max. pressure | 7 bars | min. pressure | 3 bars |
| | nebulizer pressure | 2 bars | nebulization time | 600ms | |
| | dose per nebulization | 10 ± 1 µL | with closed vent | 14 ± 1 µL | with open vent |

list of parameters

| | |
|---|--|
| Spirometry | VC, IVC, IC, ERV, IRV, TV |
| Flow / Volume | FVC (ex), FEV 6, FEV 3, FEV 1, FEV 0.75, FEV 0.7, FEV 0.5, FEV1/FVC, FEV1/IVC, FEV1/VC, FEV0,5/FVC, FEV0,5/IVC, FEV1*30, PEF, MEF 25, MEF 50, MEF 75, MEF 25-MEF75, AEX, t ex, EV, FVC (in), FIV1, FIV1/VC, PIF, MIF 25, MIF 50, MIF 75, MEF50/MIF50 |
| MVV | MVV, FMVV |
| Resistance by interruption | R int |
| BGA values entered manually | pH, pCO2, pO2, BE, HCO3, Hb, SO2, Na+,K+, Ca++ |
| In case of online data transfer from analyzer | list of parameters depending on type and brand of blood gas analyzer |
| External parameters | up to 14 additional parameters, this parameter can be user defined |
| Rhinomanometry | R75, R150, R300, L75, L150, L300, RES-L75, RES-L150, RES-L300, RES-R75, RES-R150, RES-R300 |
| Respiratory drive | P0,1, PI max, PE max, P0,1 max, VE, P0,1 / VE, P0,1 / PI max, VT, BF, VT/Ti, Ti/Ttot, P0,1/VT/Ti, PI max 1.0, PE max 1.0 |
| Compliance | C static, C dynamic |
| Resistance measured by bodyplethysmography | RAW (eff), sRAW, GAW, sGAW |
| Volumes measured by bodyplethysmography | TGV, RV, TLC, TGV/TLC, RV/TLC, IC, VC, IC/TLC, VC/TLC, IC/TGV, VC/TGV, IC/RV, VC/RV |
| CO diffusion | TLco, Kco, TLC, RV, FRC, VA, t diff, Hb, VIN (ch4) |

legal notice


We accept no liability for typographical errors or technical errors or content errors in this complete document. Pictures may not represent the current state. We assume no guarantee for the availability of the products and/or features and/or options described.

The document includes trademark of Microsoft Corporation and the Medical Equipment Europe GmbH

© 2015 by Alfred Albert

About us

Medical Equipment Europe GmbH is specialized in development, manufacturing and distribution of measuring systems for cardiopulmonary function diagnostics.



Very precise and easy to use measuring technology combined with patient acceptance are main targets of our development.

We design and produce all our devices within one production place on a surface of more than 3.500 m². We dispose a high vertical range of manufacturing.

We concern to produce our products with a minimum of environmental pollution. Transport of semifinished products we keep at a minimum.



**medical
equipment
europe**

our distribution partner