


# TECHNICAL DATA

## ENBIO<sub>s</sub>

|   |                                   |
|---|-----------------------------------|
| Class                                   | <b>B</b>                          |
| Power supply                            | <b>220-230V AC / 50/60Hz</b>      |
| Installed power                         | <b>3250 W</b>                     |
| Max. power consumption                  | <b>15A</b>                        |
| Working pressure                        | <b>2,1 Bar</b>                    |
| Max. pressure                           | <b>2,3 Bar</b>                    |
| Max. process temperature                | <b>137°C</b>                      |
| Process chamber volume                  | <b>2,7 l</b>                      |
| Programs                                | <b>134°C FAST   134°C   121°C</b> |
| Max. packages amount (circa)            | <b>8</b>                          |
| Max. weight of the load                 | <b>0,5 kg</b>                     |
| Mass                                    | <b>15 kg</b>                      |
| Level of security                       | <b>IP20</b>                       |
| Data archiving                          | <b>USB Memory stick</b>           |
| Water supply*                           | <b>external</b>                   |
| * demineralized or distilled water only |                                   |



**Enbio Group AG**  
 Eichengasse 3  
 4702 Oensingen  
 Switzerland