

Technical Data Audiometer MA 25

MA 25 SPECIFICATIONS

Test Signals	Pure tone and warble tone (pulsed and continuous)
Test Frequencies	125, 250, 500, 750, 1000, 1500, 2000, 3000, 4000, 6000, 8000 Hz
Level Steps	5 or 1 dB
Hearing Level Range AC	-10 dB HL to 100 dB HL

DEVICE GENERAL

Dimensions	W x D x H: 22.5 cm x 18 cm x 5.5 cm
Weight	980 g including standard components
Power Supply	3 AA batteries or optional power supply (100-240 VAC $\pm$ 10 %, 50/60 Hz, 0.6 A)
Display	38.1 mm x 50.8 mm, monochrome
Language	English, Dutch, French, German, Italian, Spanish

STANDARDS

Safety	IEC 60601-1: 2012 AAMI ES60601-1:2005+A2+A1 CAN/CSA-C22.2 No. 60601-1:14 Type B applied parts
EMC	IEC 60601-1-2:2014
CE0123	Class IIa according to Medical Device Regulation (EU) 2017/745
Audiometer	IEC 60645-1:2017 Type 4 ANSI/ASA S3.6-2018 Type 4

STANDARD COMPONENTS

- AC headset DD45
- 3 AA batteries
- Audiogram Pad

OPTIONAL TRANSDUCER

- AC noise-enclosure headset DD65 v2

OPTIONAL ACCESSORIES

- Patient response switch
- Carrying case
- Power supply

MA 25E VERSION

Special Tests	Automated Hughson-Westlake
Additional Accessories	• Patient response switch • USB cable
Software	MAICO Sessions PC Software
Optional Software	• OtoAccess® Database • OtoAccess® Worklist HL7 (OtoAccess® Database add-on) • Noah Database

Specifications are subject to change without notice.